

We make the tools.  
You'll win  
the fight.



## VIROTROL SARS-CoV-2 and VIROCLEAR SARS-CoV-2

Serological Controls for SARS-CoV-2 Antibody Testing

**BIO-RAD**

## VIROTROL SARS-CoV-2

### Reactive for SARS-CoV-2 total IgG/IgM and IgG antibodies

An unassayed, reactive serological quality control for the qualitative determination of total IgG/IgM and IgG antibodies to SARS-CoV-2.

- IVD/CE Mark
- Unassayed reactive control
- Human plasma based
- 2 year shelf life at 2-8°C
- 60 day open-vial stability at 2-8°C

### Analytes

SARS-CoV-2 Total IgG/IgM  
SARS-CoV-2 IgG

## VIROCLEAR SARS-CoV-2

### Non-reactive for SARS-CoV-2 total IgG/IgM and IgG antibodies

An unassayed, non-reactive serological quality control for the qualitative determination of total IgG/IgM and IgG antibodies to SARS-CoV-2.

- IVD/CE Mark
- Unassayed non-reactive control
- Human plasma based
- 2 year shelf life at 2-8°C
- 60 day open-vial stability at 2-8°C

### Analytes

SARS-CoV-2 Total IgG/IgM  
SARS-CoV-2 IgG

### List of Popular Instruments, Assays, Methods and Bio-Rad Controls

| Manufacturer/ Instrument        | Assay/Method                   | Product Name<br>Catalog Numbers | VIROTROL SARS-CoV-2 |                    |                    | VIROCLEAR SARS CoV-2 |
|---------------------------------|--------------------------------|---------------------------------|---------------------|--------------------|--------------------|----------------------|
|                                 |                                |                                 | 200300A<br>200305A  | 200300B<br>200305B | 200300C<br>200305C | 200500<br>200505     |
| Abbott ARCHITECT i System       | SARS-CoV-2 IgG                 |                                 | ●                   |                    |                    | ●                    |
| Babson Diagnostics              | aC19G1 SARS-CoV-2 IgG          |                                 | ●                   |                    |                    | ●                    |
| Beckman Coulter                 | Access SARS-CoV-2 IgG          |                                 | ●                   |                    |                    | ●                    |
| bioMérieux VIDAS                | Anti-SARS-CoV-2 IgG            |                                 | ●                   | ●                  |                    | ●                    |
| Bio-Rad Laboratories            | Platelia SARS-CoV-2 Total Ab   |                                 | ●                   |                    |                    | ●                    |
| Cellex Rapid Test               | qSARS-CoV-2 IgG/IgM Rapid Test |                                 | ●                   |                    |                    | ●                    |
| Diasorin Liaison                | SARS-CoV-2 S1/S2 IgG           |                                 | ●                   |                    |                    | ●                    |
| EuroImmun                       | Anti-SARS-CoV-2 ELISA IgG      |                                 | ●                   |                    |                    | ●                    |
| Ortho VITROS                    | Anti-SARS-CoV-2 IgG            |                                 | ●                   |                    |                    | ●                    |
| ECi/ECiQ/3600/5600/XT/7600      | Anti-SARS-CoV-2 Total          |                                 |                     |                    | ●                  | ●                    |
| Roche Cobas e411/e601/e602/e801 | Elecsys Anti-SARS-CoV-2        |                                 |                     | ●                  |                    | ●                    |
| Siemens Atellica IM             | SARS-CoV-2 Total (COV2T)       |                                 |                     | ●                  |                    | ●                    |
| Siemens ADVIA Centaur           | SARS-CoV-2 Total               |                                 |                     | ●                  |                    | ●                    |

The controls listed here include analyte(s) that may be appropriate for use on the test systems noted.

Content is provided as a guide and for informational purposes only. Results and specific performance characteristics may vary with different manufacturers' test kits and assay procedures, different lots of the same test kit, and different lots of the control material. Please note that not all products may be available in all countries. Refer to the product insert for current product claims.

## Ordering Information

### Catalog No. Description

#### VIROTROL SARS-CoV-2

|         |                             |          |
|---------|-----------------------------|----------|
| 200300A | Single Level/Tube . . . . . | 1 x 4 mL |
| 200305A | Single Level/Tube . . . . . | 5 x 4 mL |
| 200300B | Single Level/Tube . . . . . | 1 x 4 mL |
| 200305B | Single Level/Tube . . . . . | 5 x 4 mL |
| 200300C | Single Level/Tube . . . . . | 1 x 4 mL |
| 200305C | Single Level/Tube . . . . . | 5 x 4 mL |

### Catalog No. Description

#### VIROCLEAR SARS-CoV-2

|        |                             |          |
|--------|-----------------------------|----------|
| 200500 | Single Level/Tube . . . . . | 1 x 4 mL |
| 200505 | Single Level/Tube . . . . . | 5 x 4 mL |

Refer to myeinserts.com or the package insert of currently available lots for specific analyte and stability claims.



**Bio-Rad  
Laboratories, Inc.**

**For further information, please contact the Bio-Rad office nearest you or visit our website at [www.bio-rad.com/qualitycontrol](http://www.bio-rad.com/qualitycontrol)**

Clinical  
Diagnostics Group

Website [www.bio-rad.com/qualitycontrol](http://www.bio-rad.com/qualitycontrol) U.S. 1 800 224 6723 Australia +61 (2) 9914 2800 Austria +43 (0) 1 877 89 01 9 Belgium +32 (0) 3 710 53 00 Brazil +55 11 3065 7550 Canada +1 514 334 4372 China +86 21 6169 8500 Czech Republic +420 241 431 660 Denmark +45 44 52 10 00 Finland +358 9 804 22 00 France +33 (0) 1 47 95 60 00 Germany +49 (0) 89 31884 393 Greece +30 210 7774396 Hong Kong +85 2 2789 3300 Hungary +36 1 459 6190 India +91 124 4029300 Israel +972 03 963 6025 Italy +39 024 94 86 600 Japan +81 3 6361 7070 Korea +82 080 007 7373 Mexico +52 (55) 5488 7670 The Netherlands +31 (0) 318 540 666 New Zealand +64 (9) 415 2280 Norway +47 23 38 41 30 Poland +48 22 331 99 99 Portugal +351 21 47 27 700 Russia +7 495 721 1404 Singapore +65 6415 3170 South Africa +27 11 442 8508 Spain +34 91 490 6580 Sweden +46 844 98053 Switzerland +41 (0) 61 717 9555 Taiwan +886 (2) 2578 7189 Thailand (662) 651 8311 United Kingdom +44 (0) 1923 471301

BIO-RAD, VIROCLEAR and VIROTROL are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner.